

ELISA KITS & REAGENTS



ELISAs

25 different milk and serum proteins

ADAMTS13
Albumin
Alpha-1-Antitrypsin
Alpha-Lactalbumin
Beta-2 Glycoprotein 1
Beta-Lactoglobulin
CRP
Hemoglobin
Hemoglobin F
IgA
IgD
IgE
IgG
IgG(T)
IgG1
IgG2
IgG2a
IgG2b
IgG2c
IgG3
IgM
Kappa
Lactoferrin
Lambda
Transferrin

across 10 species

- Bovine
- Chicken
- Dog
- Goat
- Horse
- Human
- Mouse
- Pig
- Rabbit
- Rat

in 2 formats

ELISA Kits - Kits include pre-coated 96-well strip plate, standard, detection antibody, sample dilution buffer, HRP solution, TMB substrate, stop solution, 20X wash buffer, sealing tape and detailed instructions.

ELISA Sets - Affinity purified antibody, calibrator and HRP conjugate are packaged to provide a quantitative ELISA Set. Each set contains reagents sufficient for 1000 single well assays. Sets do not include plates, buffers or substrate. If these components are not readily available in your lab, it is recommended to consider ordering an ELISA Starter Accessory Kit, Catalog No. E101 or Catalog No. E103.

ELISA Accessories

| Description | Cat. No. | Description | Cat. No. |
|--------------------------------|----------|---------------------|----------|
| ELISA Blocking Buffer | E104 | ELISA Wash Solution | E106 |
| ELISA Coating Buffer | E107 | Microtiter Plate | E105 |
| ELISA Starter Accessory Kit | E101 | TMB Substrate | E102 |
| ELISA Starter Accessory Kit II | E103 | 10% Tween 20 | E108 |
| ELISA Stop Solution | E115 | | |




ELISA Kits

Kits include pre-coated 96-well strip plate, standard, detection antibody, sample dilution buffer, HRP solution, TMB substrate, stop solution, 20X wash buffer, sealing tape and detailed instructions.

| Description | Range | Cat. No. | |
|---------------------------------|-------------------|----------|---------|
| Bovine Albumin ELISA Kit | 0.69 – 500 ng/ml | E11-113 | BOVINE |
| Bovine IgA ELISA Kit | 1.37 – 1000 ng/ml | E11-121 | |
| Bovine IgG ELISA Kit | 0.69 – 500 ng/ml | E11-118 | |
| Bovine IgG1 ELISA Kit | 1.0 – 750 ng/ml | E11-116 | |
| Bovine IgG2 ELISA Kit | 1.0 – 750 ng/ml | E11-117 | |
| Bovine IgM ELISA Kit | 1.37 – 1000 ng/ml | E11-101 | |
| Bovine Lactoferrin ELISA Kit | 0.69 – 500 ng/ml | E11-126 | |
| Chicken IgA ELISA Kit | 1.37 – 1000 ng/ml | E33-103 | CHICKEN |
| Chicken IgG ELISA Kit | 0.69 – 500 ng/ml | E33-104 | |
| Chicken IgM ELISA Kit | 0.69 – 500 ng/ml | E33-102 | |
| Human ADAMTS13 ELISA Kit | 12.5% – 100% | E88-500 | HUMAN |
| Human Albumin ELISA Kit | 0.69 – 500 ng/ml | E88-129 | |
| Human Alpha-1-Antitrypsin Kit | 0.24 – 200 ng/ml | E88-122 | |
| Human Beta-2 Glycoprotein 1 Kit | 1.0 – 750 ng/ml | E88-142 | |
| Human CRP ELISA Kit | 0.07 – 50 ng/ml | E88-125 | |
| | | | |

| Description | Range | Cat. No. | |
|-----------------------------|-------------------|----------|-------|
| Human Hemoglobin ELISA Kit | 0.27 – 200 ng/ml | E88-135 | HUMAN |
| Human IgA ELISA Kit | 1.0 – 750 ng/ml | E88-102 | |
| Human IgE ELISA Kit | 0.69 – 500 ng/ml | E88-108 | |
| Human IgG ELISA Kit | 0.69 – 500 ng/ml | E88-104 | |
| Human IgM ELISA Kit | 1.0 – 750 ng/ml | E88-100 | |
| Human Kappa ELISA Kit | 1.37 – 1000 ng/ml | E88-115 | |
| Human Lactoferrin ELISA Kit | 1.56 – 100 ng/ml | E88-143 | |
| Human Lambda ELISA Kit | 1.37 – 1000 ng/ml | E88-116 | |
| Human Transferrin ELISA Kit | 0.69 – 500 ng/ml | E88-128 | |
| | | | |
| Mouse Albumin ELISA Kit | 1.23 – 900 ng/ml | E99-134 | MOUSE |
| Mouse IgA ELISA Kit | 1.37 – 1000 ng/ml | E99-103 | |
| Mouse IgE ELISA Kit | 1.37 – 1000 ng/ml | E99-115 | |
| Mouse IgG ELISA Kit | 0.69 – 500 ng/ml | E99-131 | |
| Mouse IgG1 ELISA Kit | 1.0 – 750 ng/ml | E99-105 | |
| Mouse IgG2a ELISA Kit | 0.69 – 500 ng/ml | E99-107 | |
| Mouse IgG2b ELISA Kit | 0.34 – 250 ng/ml | E99-109 | |
| Mouse IgG2c ELISA Kit | 0.69 – 500 ng/ml | E99-136 | |
| Mouse IgG3 ELISA Kit | 1.37 – 1000 ng/ml | E99-111 | |
| Mouse IgM ELISA Kit | 1.37 – 1000 ng/ml | E99-101 | |
| Mouse Transferrin ELISA Kit | 0.69 – 500 ng/ml | E99-129 | |
| Pig Albumin ELISA Kit | 1.23 – 900 ng/ml | E101-110 | PIG |
| Pig IgA ELISA Kit | 1.37 – 1000 ng/ml | E101-102 | |
| Pig IgG ELISA Kit | 1.37 – 1000 ng/ml | E101-104 | |
| Pig IgM ELISA Kit | 1.37 – 1000 ng/ml | E101-100 | RBT |
| Rabbit IgG ELISA Kit | 0.69 – 500 ng/ml | E121-111 | |
| | | | |
| Rat Albumin ELISA Kit | 1.37 – 1000 ng/ml | E111-125 | RAT |
| Rat IgA ELISA Kit | 1.37 – 1000 ng/ml | E111-102 | |
| Rat IgG ELISA Kit | 1.37 – 1000 ng/ml | E111-128 | |
| Rat IgM ELISA Kit | 2.74 – 2000 ng/ml | E111-100 | |



ELISA Sets

Affinity purified antibody, calibrator and HRP conjugate are packaged to provide a quantitative ELISA Set. Each set contains reagents sufficient for 1000 single well assays. Sets do not include plates, buffers or substrate. If these components are not readily available, it is recommended to consider ordering an ELISA Starter Accessory Kit, Catalog No. E101 or Catalog No. E103.

| Description | Cat. No. | |
|--|----------|---------|
| Bovine Albumin ELISA Quantitation Set | E10-113 | BOVINE |
| Bovine alpha-Lactalbumin ELISA Quantitation Set | E10-128 | |
| Bovine beta-Lactoglobulin ELISA Quantitation Set | E10-125 | |
| Bovine IgA ELISA Quantitation Set | E10-121 | |
| Bovine IgG1 ELISA Quantitation Set | E10-116 | |
| Bovine IgG2 ELISA Quantitation Set | E10-117 | |
| Bovine IgG ELISA Quantitation Set | E10-118 | |
| Bovine IgM ELISA Quantitation Set | E10-101 | |
| Bovine Lactoferrin ELISA Quantitation Set | E10-126 | |
| Bovine Transferrin ELISA Quantitation Set | E10-122 | |
| Chicken IgA ELISA Quantitation Set | E30-103 | CHICKEN |
| Chicken IgG ELISA Quantitation Set | E30-104 | |
| Chicken IgM ELISA Quantitation Set | E30-102 | |
| Dog IgA ELISA Quantitation Set | E40-104 | DOG |
| Dog IgE ELISA Quantitation Set | E40-125 | |
| Dog IgG ELISA Quantitation Set | E40-118 | |
| Dog IgM ELISA Quantitation Set | E40-116 | GOAT |
| Goat IgA ELISA Quantitation Set | E50-106 | |
| Goat IgG ELISA Quantitation Set | E50-104 | |
| Goat IgM ELISA Quantitation Set | E50-105 | HORSE |
| Horse IgA ELISA Quantitation Set | E70-116 | |
| Horse IgG(T) ELISA Quantitation Set | E70-105 | |
| Horse IgM ELISA Quantitation Set | E70-114 | HUMAN |
| Human Albumin ELISA Quantitation Set | E80-129 | |
| Human Hemoglobin ELISA Quantitation Set | E80-135 | |
| Human Hemoglobin F ELISA Quantitation Set | E80-136 | |
| Human IgA ELISA Quantitation Set | E80-102 | |
| Human IgD ELISA Quantitation Set | E80-106 | |
| Human IgE ELISA Quantitation Set | E80-108 | |

| Description | Cat. No. | |
|--|----------|-------|
| Human IgG ELISA Quantitation Set | E80-104 | HUMAN |
| Human IgM ELISA Quantitation Set | E80-100 | |
| Human Kappa ELISA Quantitation Set | E80-115 | |
| Human Lambda ELISA Quantitation Set | E80-116 | |
| Human Transferrin ELISA Quantitation Set | E80-128 | |
| Mouse Albumin ELISA Quantitation Set | E90-134 | |
| Mouse IgA ELISA Quantitation Set | E90-103 | |
| Mouse IgE ELISA Quantitation Set | E90-115 | |
| Mouse IgG ELISA Quantitation Set | E90-131 | |
| Mouse IgG1 ELISA Quantitation Set | E90-105 | |
| Mouse IgG2a ELISA Quantitation Set | E90-107 | |
| Mouse IgG2b ELISA Quantitation Set | E90-109 | |
| Mouse IgG2c ELISA Quantitation Set | E90-136 | |
| Mouse IgG3 ELISA Quantitation Set | E90-111 | |
| Mouse IgM ELISA Quantitation Set | E90-101 | PIG |
| Pig Albumin ELISA Quantitation Set | E100-110 | |
| Pig IgA ELISA Quantitation Set | E100-102 | |
| Pig IgG ELISA Quantitation Set | E100-104 | RBT |
| Pig IgM ELISA Quantitation Set | E100-100 | |
| Rabbit IgG ELISA Quantitation Set | E120-111 | RAT |
| Rabbit IgM ELISA Quantitation Set | E120-110 | |
| Rat Albumin ELISA Quantitation Set | E110-125 | |
| Rat IgA ELISA Quantitation Set | E110-102 | |
| Rat IgG ELISA Quantitation Set | E110-128 | |
| Rat IgG1 ELISA Quantitation Set | E110-106 | |
| Rat IgG2a ELISA Quantitation Set | E110-109 | |
| Rat IgG2b ELISA Quantitation Set | E110-111 | |
| Rat IgG2c ELISA Quantitation Set | E110-114 | |
| Rat IgM ELISA Quantitation Set | E110-100 | |



For *in vitro* laboratory use only. Not for any clinical, therapeutic or diagnostic use in humans or animals. Not for animal or human consumption.

Warranty : Products are warranted by Bethyl Laboratories, Inc. to meet stated product specifications and to conform to label descriptions when used, handled and stored according to instructions. Unless otherwise stated, this warranty is limited to six months from date of sale. Bethyl Laboratories sole liability for the product is limited to replacement of the product or refund of the purchase price. Bethyl Laboratories products are supplied for research applications. They are not intended for medicinal, diagnostic or therapeutic use. The products may not be resold, modified for resale or used to manufacture commercial products without prior written approval from Bethyl Laboratories, Inc.